


Signing Off a Result

Once a Provider has reviewed a result, it can be signed off and removed from the Investigation Viewer.

To Sign Off a Result:

1. Access the Investigation Viewer and make sure the **Results** tab is displayed.
2. Click on the Lab Report you want to sign-off.
3. Click the  icon.

The result is removed from the Investigation Viewer.

Notice the **Signed Off by** box at the bottom of the Investigation Viewer. It shows the date and time the result was signed-off, the Provider's ID and the Provider's ABELMed log on User Name.

Example:

Signed off by:			
24/Nov/11	19:44	PC	Dr. PaulChiro-554321